

Modelo 240

Audiómetro de diagnóstico portátil y automático

El Amplivox 240 es un audiómetro de diagnóstico rentable, que es ideal tanto para consultorio y el uso móvil.

Funciones

El Amplivox 240 ofrece la conducción aérea y ósea. Las características estándar incluyen una función automática de prueba (AC), el enmascaramiento de banda estrecha, un rango de atenuador de hasta 120 dB HL y una función talkover integral.

Facil de usar

El panel de control es intuitivo hace las pruebas de diagnóstico lo más simple y directo posible.

La función de retención de umbral reduce el tiempo de prueba y da un mayor nivel de concentración sujeto de prueba.

- Fácil de utilizar
- Compacto y portátil
- Conectividad USB
- Función de retención de Umbral
- Auriculares de inserción opcionales
- 11 frecuencias de prueba
- -10 a 120 dBHL Vía aérea
- Prueba aérea automática

Totalmente portátil

Diseñado para satisfacer todos los requisitos de pruebas domiciliarias, el Amplivox 240 cuenta con adaptador de alimentación de red estándar. Su tamaño compacto, peso mínimo y opciones de energía hacen de este un audiómetro de diagnóstico totalmente portátil.



especialistas en productos para la salud

Model 240

Portable diagnostic audiometer



Technical specifications

Output/input data

| | |
|----------------------------|---|
| Outputs: | Left AC, Right AC, Bone (L&R) insert masking |
| Frequency range (Hz): | AC: 125-8kHz BC: 250Hz-8kHz |
| Frequency accuracy: | <1% |
| Output level range (AC): | -10dBHL to 120dBHL (max) |
| Output level range (BC): | -10dBHL to 70dBHL (max) |
| Insert masking output: | 90dBHL max (250-4kHz) |
| Output level accuracy: | Within 3dB |
| Output level step size: | 5dB |
| Output transducer options: | Standard audiometric headset and/or optional insert earphones |
| Tone present: | Single, pulsed or warble |
| AC test: | Manual or automatic mode |
| Masking: | Narrowband |
| Communication: | Integral talk forward facility |

Physical data

| | |
|------------------|------------------------------------|
| Display: | 2 lines of 24 characters (backlit) |
| Power (battery): | 4 x 1.5V "C" type batteries |
| Power (mains): | 100-240V ac 50-60Hz |
| Dimensions (mm): | 270 x 175 x 68 high |
| Weight: | 735 grams |

Standards

The Amplivox 240 conforms to the relevant clauses of the following standards.

| | |
|--------------|--|
| Safety: | IEC 60601-1 |
| EMC: | IEC 60601-1-2 |
| Audiometric: | IEC 60645-1 |
| CE mark: | Complies to EU Medical Device Directive |

Equipamiento

- Auriculares
- Auriculares de conducción ósea
- Interruptor de respuesta del paciente
- Tarjetas de audiometría (50)
- Manual de instrucción
- Maletín de transporte



Carrying case

For domiciliary use, the Amplivox 240 is supplied with a carrying case, designed to accommodate the audiometer and all accessories.

Optional equipment

- Masking earpiece and lead
- Audiocups (noise reducing headset enclosures)
- Insert earphones

amplivox

Amplivox Ltd, 29-30 Station Approach
Kidlington, Oxford, United Kingdom OX5 1JD
Tel: +44 (0)1865 842411 • Fax: +44 (0)1865 841853
E-mail: sales@amplivox.ltd.uk

www.amplivox.ltd.uk